

One Platform to Power the Future of Aesthetic Treatments.

The Icon™ incorporates laser and Optimized Pulsed Light into one platform to deliver a comprehensive suite of treatments that patients want.





Facial Skin Revitalization – MaxG™ OPL™

- Optimized for clearing vascular & pigmented lesions
- Dynamic Spectrum Shifting[™] distributes energy across highly selective
 spectra to optimize closure of all types of vessels including larger deeper facial vessels
- Unique dual-band filters optimize absorption in hemoglobin and epidermal melanin while minimizing heating, providing a comfortable no-downtime treatment





Photos courtesy of M. Sinclair, MD

Photos courtesy of H. Alibhai, MD

Fractional Non-Ablative Skin Resurfacing – **1540 Laser**

- Improves skin irregularities and reduces the appearance of acne scars and surgical scars
- Clinically proven laser for the treatment of stretch marks
- Customizable facial treatments to focus on each patient's problem area



XF Microlens™ high speed with increased surface coverage



XD Microlens™

uses point compression and a unique lens design to allow for greater depth of penetration



Photos courtesy of S. Doherty, MD




Photos courtesy of D. Vasily, MD

Pigmented Lesion Removal – **MaxYs[™] OPL[™]**

- Delivers optimally filtered light for the treatment of unwanted pigmented lesions; ideal for chest, arms and hands
- · High peak power to remove difficult to treat finer, lighter hair









Photos courtesy of Cynosure, Inc.



Permanent Hair Reduction – MaxR™ & MaxRs™ OPL™

- Hair removal on all skin types including tanned skin
 - o MaxR faster coverage for large areas including back, legs, and chest
 - o MaxRs perfect for small areas including face, underarms, and bikini





Photos courtesy of K. Khatri, MD

Photos courtesy of K. Khatri, MD



Fractional Ablative Skin Resurfacing – 2940 Laser

- Perfect solution for textural irregularities and wrinkles
- One-time treatment that provides dramatic results
- Less discomfort and downtime than traditional CO₂ lasers
- Customizable treatments to meet desired results and downtime





Photos courtesy of D. Kent, MD

Photos courtesy of H. Alibhai, MD • 1540/2940 combination



Leg Vein Removal - 1064+ Laser

- Comfortably clears a variety of leg vessels including smaller red vessels and larger blue reticular leg veins
- · Additional benefits include permanent reduction of unwanted hair on skin types V and VI





Photos courtesy of D. Vasily, MD

Photos courtesy of D. Vasily, MD



Handpiece	Spot Size (mm)	Spectral Rang	je (nm)	Fluence Range* (J/cm²)
MaxG	10 x 15	500-670 & 870	-1200	Up to 80
MaxYs	12 x 28	525-1200		Up to 90
MaxR	16 x 46	650-1200		Up to 48
MaxRs	12 x 28	650-1200		Up to 72
Handpiece	Wavelength (nm)	Laser Type	Spot Size (mm)	Pulsewidth (ms)
Handpiece 1540	Wavelength (nm) 1540	Laser Type Er:Glass	Spot Size (mm) 12 & 15	Pulsewidth (ms) 10 & 15
_		21		` '

^{*}The effective fluence may be higher due to photon recycling, wide spectrum, spot size and skin contact.



Corporate Headquarters

5 Carlisle Road Westford, Massachusetts 01886 USA

TELEPHONE 978 256 4200 TOLL FREE 800 886 2966 www.cynosure.com







© 2015 Cynosure, Inc. Cynosure and Palomar are registered trademarks of Cynosure, Inc. Icon, MaxG, XF Microlens, XD Microlens, MaxYs, MaxRs, are trademarks of Cynosure, Inc. 921-0532-000 Rev. 1 04/15